GENERAL

The Powerex laboratory vacuum system is designed to create a suction system to remove unwanted fluids or gases from laboratory working areas. The laboratory vacuum system package is compliant with the NFPA 99 requirements for Category 3 laboratory systems. Each system is completely tested before shipment and includes:

- Multiple vacuum pumps and associated equipment
- ASME air receiver
- NEMA 12 multiplex control panel

Each pump is factory piped to a common intake manifold. Vibration isolation pads are included with the system.

ROTARY VANE VACUUM PUMP

Each vacuum pump shall be of the rotary vane aircooled design. Each vacuum pump shall be directdriven through a shaft coupling by a TEFC electric motor. Belt drives shall not be permitted. The vacuum pumps shall be mounted on vibration isolators.

- Each vacuum pump shall be air-cooled and have absolutely no water requirements.
- Each pump is equipped with a 5 micron inlet filter for removal of particulates.
- Lubrication shall be provided by an integral, fully recirculating oil supply that is filtered by an automotive-type, spin-on oil filter. Non-recirculating (once-through) or partial re-circulating oil supply systems shall not be permitted. Each vacuum pump is supplied with ISO100 PAO synthetic oil.
- The oil separation systems shall be integral and shall consist of no less than three stages of internally installed oil and smoke eliminators through which the exhaust gas stream must pass. This system shall consist of bulk separation, oil mist elimination, and smoke elimination, and shall be capable of removing 99.9+ percent of all oil and smoke particles from the exhaust gas stream.
- Each vacuum pump shall include a built-in, antisuck-back valve, mounted at the pump inlet, and three sliding vanes.

Each pump is initially filled with full synthetic oil for improved lubrication, longer life, and chemical resistance. Laboratory Lubricated Vane Vacuum System with Premium NFPA Controls Rev. 02/27/14

<u>MOTOR</u>

The motor is continuous duty, C-face, TEFC, suitable for 208-230, or 460V, 3 phase, 60 hertz electrical operation.

AIR RECEIVER

The system shall include an ASME rated air receiver. The tank shall be internally lined and come with a manual drain and sight gauge. Systems with Basic Controls are also equipped with a vacuum gauge. Receiver is plumbed so that system vacuum must pass through the air receiver prior to the pumps to provide additional protection for the pumps.

PREMIUM NFPA CONTROL PANEL

The control system provides automatic lead/lag sequencing and automatic alternation of all pumps in order to equalize the amount of usage among the available vacuum pumps. The Premium NFPA Control Panel shall include a gateway server card and all features listed below:

- PLC controller and a 6" color touch screen panel which displays the operating status of the unit
- Building automation communication gateway with BacNet[®] protocol and Web server features. Web server features include email notifications in case the system is in alarm or has achieved one its maintenance intervals and requires service.
- Ethernet port for connection to BacNet[®] server or direct connection to facility Ethernet for viewing of system operations and status via device IP-address
- UL508A listed control panel in a NEMA 12 enclosure. The panel door will include: the color touch screen, audible and visual alarms with an acknowledge button, and an HOA switch for each pump
- Magnetic starters
- Vacuum transducer for process control
- Single point power connection
- Redundant 120Vac control transformers with fused primary and secondary protection
- System overload trip, high temperature conditions, low vacuum or maintenance intervals for the pump will result in visual and audible alarms.

OPTIONAL AUTO PURGE

NOTE: This option is available with our Premium NFPA Control Panel. It is not available on systems with Basic NFPA controls.

Laboratory Lubricated Vane Vacuum System with Premium NFPA Controls Rev. 02/27/14

Systems with the auto purge feature will include:

- Control and operation of air actuated and/or electric isolation and purge valves
- 3-position mode switch (Auto-Off-Manual)
- Pump Status (PUMP ON, PUMP READY, PUMP OFF, PUMP PURGING)

Pumps will run through an automatic purge cycle every 24 hours (after alternation phase), to purge any gases from the pump to atmosphere. A 3-position switch is provided to select automatic purge, manual purge, or purge off modes. (**NOTE**: Pneumatic purge isolation valves require 80 psi compressed air to be field supplied.)




MODEL

DIMENSIONS

Laboratory Vertical **Duplex Lubricated Vane** Vacuum Package with **Premium Controls**

Rev. 01/08/13




	Laboratory Lubricated Vane Vacuum Package with Premium Controls												
		Tank Size	SCFM @		dB(A)	SYSTEM F.L.A. ²			SYSTEM				
MODEL	HP	(Gal.)	19" Hg ²	BTU/HR ²	LEVEL	208V	230V	460V	WEIGHT (LBS.)				
LVVD0203	2 (2)	80 V	22	8,670	73	13.2	12.8	7.4	800				
LVVD0303	3 (2)	80 V	34	13,006	73	18.8	17.6	9.8	1,000				

Notes:

1-3 Year Limited Warranty



Duplex Horizontal Laboratory Lubricated Vane Vacuum Package with Premium Controls



	Laboratory Lubricated Vane Vacuum Package with Premium Controls												
		Tank Size	SCFM @		dB(A)	SYSTEM F.L.A. ²		SYSTEM WEIGHT					
MODEL	HP	(Gal.)	19" Hg ²	BTU/HR ²	LEVEL	208V	230V	460V	(LBS.)				
LVD0203	2 (2)	80 H	22	8,670	73	13.2	12.8	7.4	950				
LVD0303	3 (2)	80 H	34	13,006	73	18.8	17.6	9.8	1,100				
LVD0404	5 (2)	120 H	52	21,676	74	29.6	28.0	15.0	1,250				
LV D0504	5 (2)	120 H	82	21,676	77	30.4	28.0	15.0	1,250				


1-3 Year Limited Warranty

Notes:




Duplex Laboratory Lubricated Vane Vacuum Package with Premium Controls

DIMENSIONS DIM. DIM. C MODEL DIM. A в Outlet Inlet 41" 68' 95" LVPDT0402 2" 2" LVPDT0502 39' 68" 95" 2" 2" LVPDT0752 42' 74" 95' 3" 3"



·R



	Laboratory Lubricated Vane Vacuum Package with Premium Controls												
		Tank Size	SCFM @		dB(A)	SYSTEM F.L.A. ²			SYSTEM WEIGHT				
MODEL	HP	(Gal.)	19" Hg ²	BTU/HR ²	LEVEL	208V	230V	460V	(LBS.)				
LVPDT0402	5 (2)	60 H*	52	21676	74	29.6	28.0	15.0	1,600				
LVPDT0502	5 (2)	60 H*	74	21,676	77	30.0	29.6	15.8	1,815				
LVPDT0752	7.5 (2)	60 H*	104	32,512	82	44.0	42.8	22.4	2,300				

1-3 Year Limited Warranty

Notes:


2 - SCFM, BTU and FLA values listed are with all pumps running

** Receiver is frame mounted on top of the two pumps.




Duplex Laboratory Lubricated Vane Vacuum Package with **Premium Controls**

Rev. 03/22/13



			DIM ENSI	ONS		
	MODEL	DIM. A	DIM. B	DIM. C	Inlet	Outlet (x2)
ĺ	LVPD1005	66"	59"	83"	3"	2"
ĺ	LVPD1505	90"	75"	88"	4"	3"
	LV PD2005	90"	75"	88"	4"	3"
	LV PD2505	90"	75"	88"	4"	3"



B

	Laboratory Lubricated Vane Vacuum Package with Premium Controls												
		Tank Size	SCFM @		dB(A)	SYSTEM F.L.A. ²			SYSTEM				
MODEL	HP	(Gal.)	19" Hg ²	BTU/HR ²	LEVEL	208V	230V	460V	WEIGHT (LBS.)				
LVPD1005	10 (2)	200 V	154	43,350	84	54.8	53.2	27.6	2,600				
LVPD1505	15 (2)	200 V	222	65,026	86	82.8	79.6	40.8	3,200				
LVPD2005	20 (2)	200 V	274	86,700	87	110.0	106.0	54.0	3,600				
LVPD2505	25 (2)	200 V	336	108,376	88	138.0	126.0	64.0	5,000				

Notes:

1-3 Year Limited Warranty



MODEL

0 0


٥

0

DIMENSIONS

Triplex Laboratory Lubricated Vane Vacuum Package 5-10 HP with **Premium Controls**

Rev. 03/22/13



	Laboratory Lubricated Vane Vacuum Package with Premium Controls												
		Tank Size	SCFM @		dB(A)	SYSTEM F.L.A. ²			SYSTEM WEIGHT				
MODEL	HP	(Gal.)	19" Hg ²	BTU/HR ²	LEVEL	208V	230V	460V	(LBS.)				
LVPT0505	5 (3)	200 V	111	32,514	79	44.0	43.4	22.7	2,100				
LVPT0755	7.5 (3)	200 V	156	48,768	84	65.0	63.2	32.6	2,750				
LVPT1005	10 (3)	200 V	231	65,025	86	81.2	78.8	40.4	3,050				

Notes:

1-3 Year Limited Warranty 2 - SCFM, BTU and FLA values listed are with all pumps running

Powerex - 150 Production Dr. - Harrison, OH 45030 - Phone: 888-769-7979 - Fax: 513-367-3125

Outlet

(X3)

3"

3"

3"

Inlet

4"

4"

4"



DIM. A

135'

135

135'

MODEL

LVPT1505 LVPT2005

LVPT2505

DIMENSIONS

88"

88"

88'

DIM. BIDIM.C


75"


75'


75"

Triplex Laboratory Lubricated Vane Vacuum Package 15-25 HP with Premium Controls

Rev. 01/08/13







	Laboratory Lubricated Vane Vacuum Package with Premium Controls												
		Tank Size	SCFM @		dB(A)	SYSTEM F.L.A. ²			SYSTEM WEIGHT				
MODEL	HP	(Gal.)	19" Hg ²	BTU/HR ²	LEVEL	208V	230V	460V	(LBS.)				
LVPT1505	15 (3)	200 V	333	97,539	88	123.2	118.4	60.2	5,800				
LVPT2005	20 (3)	200 V	411	130,050	89	164.2	158.0	80.0	6,300				
LVPT2505	25 (3)	200 V	504	162,564	90	206.0	188.0	95.0	6,950				

Notes:

¹⁻³ Year Limited Warranty

Outlet

(X4)

2"

2"

2"

Inlet

3"

3"

4"



DIM. A

98"

98'

98"

MODEL

LVPQ0505

LVPQ0755

LVPQ1005

DIMENSIONS

DIM. BIDIM.C

59"

59'

66"


83"


83"


83"

Quadplex Laboratory Lubricated Vane Vacuum Package with Premium Controls

Rev. 03/22/13







	Laboratory Claw Vacuum Package with Premium Controls												
		Tank Size	SCFM @		dB(A)	SYS	TEM F.L	SYSTEM WEIGHT					
MODEL	HP	(Gal.)	19" Hg ²	BTU/HR ²	LEVEL	208V	230V	460V	(LBS.)				
LVPQ0505	5 (4)	200 V	148	43,352	80	58.0	57.2	29.6	2,500				
LVPQ0755	7.5 (4)	200 V	208	65,024	85	86.0	83.6	42.8	3,000				
LVPQ1005	10 (4)	200 V	308	86,700	87	107.6	104.4	53.2	3,600				


Notes:


1-3 Year Limited Warranty




Quadplex Laboratory Lubricated Vane Vacuum Package with Premium Controls

DIMENSIONS Outlet MODEL DIM. A DIM. В DIM. C Inlet (X4) LVPQ1505 135" 75" 88" 4" 3" LVPQ2005 135" 75" 88" 4" 3" LVPQ2505 135" 75" 88" 4" 3"







	Laboratory Claw Vacuum Package with Premium Controls												
		Tank Size	SCFM @		dB(A)	SYS	tem F.L	SYSTEM WEIGHT					
MODEL	HP	(Gal.)	19" Hg ²	BTU/HR ²	LEVEL	208V	230V	460V	(LBS.)				
LVPQ1505	15 (4)	200 V	444	130,052	89	163.6	157.2	79.6	7,000				
LV PQ2005	20 (4)	200 V	548	173,400	90	218.0	210.0	106.0	7,700				
LV PQ2505	25 (4)	200 V	672	216,752	91	274.0	250.0	126.0	8,600				

Notes:

1-3 Year Limited Warranty